

The Garner HDH-4000 accessory allows quick, easy and complete erasure of hard drives.



Who can benefit from the HDH-4000:

- *Police Departments*
- *Surveillance Companies*
- *IT Departments*
- *School Districts*
- *Financial Institutions*
- *Libraries*
- *Hospitals*
- *Insurance Companies*
- *and many others...*

Federal Mandates require complete erasure of Hard Drives and Tapes before Re-use or Disposal. Degaussing will allow you to meet these mandates.

HIPAA

GLB Act

FERPA

HDH-4000 Hard Drive Holder

The Garner HDH-4000 is our newest accessory. The HDH-4000 allows quick, easy and complete erasure of hard drives when used in conjunction with a Garner Eliminator CDS-4000FS conveyor belt degausser or a Garner HD-1 manual model degausser.

When used with the Eliminator CDS-4000FS the center platter rotates automatically while passing through the magnetic field. This rotation is necessary in order to completely demagnetize the high coercivity magnetic particles which store the information on the hard drive platters.

When used with the Eliminator CDS-4000FS, simply place the hard drive in the HDH-4000 and set the HDH-4000 on the Eliminator CDS-4000FS conveyor belt.

When used with the HD-1, simply set the HDH-4000 in the middle of the Erase Zone and rotate the center platter 3 full revolutions, flip the hard drive over and rotate the platter an additional 3 revolutions.

Complete Erasure Everytime...

- Place your hard drive in the HDH-4000.
- Set the HDH-4000 on the belt of the CDS-4000FS and press Erase.
- The HDH-4000 will automatically pass through the erasing field while rotating in the magnetic field and be completely erased in under 10 seconds.



HDH-4000 SPECIFICATIONS

For use with Garner Degausser Models:	Eliminator CDS-4000FS Degausser, HD-1 Degausser
Dimensions:	12.75 x 10 x 1.5 inches (326 x 256 x 38 mm)
Weight:	6 lbs (13.2 kg)